

General Cell Culturing

- Culturing cells for general use (10 cm plate)
 1. Warm cell media and trypsin
 2. Aspirate cell media
 3. Wash cells with 5 mL of PBS 3 times
 4. Add 3 mL of trypsin to plate
 5. Place cells in incubator until cells are rounded and/or floating
 6. Use a P1000 set to 1000 μ L to wash cells down off plate while it is tilted. Rotate plate to ensure all cells are washed off the plate
 7. Grab x volume of the trypsinized cells that will be kept for next generation.
 8. If the same plate is being reused:
 - (a) Keep trypsinized cells in the pipet and set aside as not to contaminate the cells or allow the tip to touch anything
 - (b) Aspirate out media from 10 cm plate
 - (c) Add cells from pipet to 10 cm plate
 - (d) Add 10 mL of appropriate cell media
 9. If a new plate is being used:
 - (a) Transfer cells to new plate
 - (b) Add 10 mL of appropriate cell media
 - (c) Throw out the old cell plate
 10. Swirl plate in a infinity symbol pattern
 11. Update plate with new information (name, date, generation, cell type) if needed
- Splitting cells for live-cell imaging (35 mm plate)
 1. Aspirate media out of plate
 2. Wash cells with 5 mL of PBS 3 times
 3. Add 3 mL of trypsin
 4. Place cells in incubator until cells are rounded
 5. Use a P1000, set to 1000 μ L , to wash cells down from tilted plate. Rotate plate to ensure all cells are washed off the plate
 6. Grab a 15 mL centrifuge tube
 7. Add 1 mL of trypsinized cells to this tube
 8. Add 1 mL of the appropriate cell media to this tube as well
 9. Mix by pipetting up and down
 10. Add 2 mL of the appropriate media to the 35 mm plate
 11. Add 40 μ L of cells to the plate, making sure to disperse the cells throughout the plate when transferring from pipet
 12. Swirl plate in a infinity symbol pattern
 13. Update plate with new information (name, date, generation, cell type)