Changing Gas Tank

- 1. Make sure all nozzles are closed on the gas regulator
- 2. Unscrew regulator from gas tank
- 3. Screw gas tank cap on top of the gas tank
 - Make sure you are orthoganol to the gas output just in case something happens
- 4. Go to gas room
- 5. Remove the chain securing the gas tank
- 6. Roll gas tank cart to the appropriate lane
- 7. Stand it upright so the plate that goes under the gas tank is on the floor in the lane of the gas tank
- 8. Carefully spin the gas tank to move it onto the gas tank cart
- 9. Secure the gas tank to the cart using the chain
- 10. Carefully lean the gas tank cart back
- 11. Transport the gas cart to its destination
- 12. Carefully lean the gas tank cart down
- 13. Remove the chain securing the gas tank
- 14. Carefully spin the gas tank to remove it from the gas tank cart
- 15. Move full gas tank to a place where it will not interfere with loading the empty gas tank
- 16. Unchain the empty gas tank from the wall
- 17. Spin the empty gas tank move it close to the gas cart
- 18. Lean the gas tank cart down
- 19. Spin the empty gas tank to get it on the gas tank cart
- 20. Secure the emtpy gas tank to the gas tank cart using the chain
- 21. Bring this emtpy gas tank to the gas room
- 22. Position the gas tank cart in front of the apppropriate lane
- 23. Lean gas tank cart down
- 24. Unchain the emtpy gas tank from the gas tank cart
- 25. Spin the empty gas tank to remove it from the gas tank cart
- 26. Place it as far back as possible (two gas tanks per lane) and chain it
- 27. Label a piece of tape with "EMPTY" and place it on the tank
- 28. Go back to location of full gas tank
- 29. Carefully spin full gas tank to its long term position
- 30. Chain the full gas tank to the wall
- 31. Unscrew gas tank cap from full gas tank and place in safe place

- 32. Remove any plasic from the nozzle
- 33. Screw the regulator initially with your hand so a wrench can be used without the regulator needing support
 - Make sure that the regulator does not interact with the nozzle on top of the gas tank
 - Make sure you are orthoganol to the gas output just in case something happens
- 34. Use a wrench to finish screwing in the regulator. Use your body weight to ensure it is as tight as possible
 - Make sure that the regulator does not interact with the nozzle on top of the gas tank
 - Make sure you are orthoganol to the gas output just in case something happens